

Product datasheet for **HP213729**

AFAF (EQTN) Human qPCR Primer Pair (NM_020641)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54586
Forward Sequence:	GGAACAGCAGTGGTCATGGATG
Reverse Sequence:	GGTCATCAACGAGATTCCCAGC
ACCN:	BC014307 , NM_020641 , NM_020641.1 , NM_020641.2 , BC014307.1 , BC144602 , BC137351 , BC137352 , BX089774
UniProt ID:	Q9NQ60
Synonyms:	AFAF; C9orf11; SPACA8
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)