

Product datasheet for **HP213746**

MRS2 Human qPCR Primer Pair (NM_020662)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57380
Forward Sequence:	CTACCGATTGGCTGACGATCTC
Reverse Sequence:	GAGCGAAAGAGAGAAGGTTCCC
ACCN:	<u>NM_020662</u> , <u>NM_020662.1</u> , <u>NM_020662.2</u> , <u>NM_020662.3</u> , <u>BC001028</u> , <u>BC018833</u> , <u>BC069009</u> , <u>BM680865</u> , <u>BM829320</u> , <u>BP338269</u>
UniProt ID:	<u>Q9HD23</u>
Synonyms:	HPT; MRS2L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)