

Product datasheet for **HP213749**

CLTRN Human qPCR Primer Pair (NM_020665)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57393
Forward Sequence:	ACCTCTCAAAGCGATGGTAGCT
Reverse Sequence:	CCACAAACCAGAATGATACCCTC
ACCN:	NM_020665 , NM_020665.1 , NM_020665.2 , NM_020665.3 , NM_020665.4 , NM_020665.5 , BC014317 , BC014317.1 , BC015099 , BC050606 , NM_020665.6
UniProt ID:	Q9HBJ8
Synonyms:	NX-17; NX17; TMEM27
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)