

Product datasheet for **HP213759**

MIF4GD Human qPCR Primer Pair (NM_020679)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57409
Forward Sequence:	ATTGCTCTCGCCTGATGCCTTC
Reverse Sequence:	CCTGCTTCCTTGCTGAACACAC
ACCN:	<u>NM_020679</u> , <u>NM_020679.1</u> , <u>NM_020679.2</u> , <u>NM_020679.3</u> , <u>BC033759</u> , <u>BC033759.1</u> , <u>BI752882</u> , <u>BU190758</u> , <u>BX404175</u>
UniProt ID:	<u>A9UHW6</u>
Synonyms:	AD023; MIFD; SLIP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)