

Product datasheet for **HP213805**

MAVS Human qPCR Primer Pair (NM_020746)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57506
Forward Sequence:	ATGGTGCTCACCAAGGTGTCTG
Reverse Sequence:	TTCAGAGCTGCTGTCTAGCCA
ACCN:	<u>NM_020746</u> , <u>NM_020746.1</u> , <u>NM_020746.2</u> , <u>NM_020746.3</u> , <u>NM_020746.4</u> , <u>BC044952</u> , <u>BC044952.1</u> , <u>BC000073</u> , <u>BC040719</u> , <u>BC040930</u> , <u>BC052583</u> , <u>BC062310</u> , <u>BC063593</u> , <u>BM834250</u> , <u>BX649084</u> , <u>NM_020746.5</u>
UniProt ID:	<u>Q7Z434</u>
Synonyms:	CARDIF; IPS-1; IPS1; VISA
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)