

## Product datasheet for **HP213834**

### PDP2 Human qPCR Primer Pair (NM\_020786)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57546
Forward Sequence:	GACATCTCGCTGGAAATCCAGG
Reverse Sequence:	CTCTTGACACCAAGGATGGCT
ACCN:	<u><a href="#">NM_020786</a></u> , <u><a href="#">NM_020786.1</a></u> , <u><a href="#">NM_020786.2</a></u> , <u><a href="#">NM_020786.3</a></u> , <u><a href="#">BC028030</a></u> , <u><a href="#">BM668681</a></u> , <u><a href="#">BU735208</a></u> , <u><a href="#">NM_020786.4</a></u>
UniProt ID:	<u><a href="#">Q9P2J9</a></u>
Synonyms:	PDPC 2; PPM2B; PPM2C2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)