

## Product datasheet for **HP213847**

### KLHL8 Human qPCR Primer Pair (NM\_020803)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57563
Forward Sequence:	GTGATGTCACACTCAAGGTTGGC
Reverse Sequence:	GCGTTTGCTTGGCTTCAGCCAT
ACCN:	<a href="#">NM_020803</a> , <a href="#">NM_020803.1</a> , <a href="#">NM_020803.2</a> , <a href="#">NM_020803.3</a> , <a href="#">NM_020803.4</a> , <a href="#">BC041384</a> , <a href="#">BC041384.2</a> , <a href="#">BC041370</a> , <a href="#">BC041901</a> , <a href="#">BX538173</a> , <a href="#">BX640727</a> , <a href="#">BX640744</a> , <a href="#">BX647454</a> , <a href="#">BX648367</a>
UniProt ID:	<a href="#">Q9P2G9</a>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)