

## Product datasheet for **HP213855**

### Dedicator of cytokinesis protein 6 (DOCK6) Human qPCR Primer Pair (NM\_020812)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57572
Forward Sequence:	TTCCTGGACGACATCACTGCCT
Reverse Sequence:	GTCACTGAGGAAGAAAGCCAGG
ACCN:	<a href="#">NM_020812</a> , <a href="#">NM_020812.1</a> , <a href="#">NM_020812.2</a> , <a href="#">NM_020812.3</a> , <a href="#">BC008335</a> , <a href="#">BC015424</a> , <a href="#">BC051330</a> , <a href="#">BC095483</a> , <a href="#">BC146786</a> , <a href="#">BI820302</a> , <a href="#">BI822534</a> , <a href="#">NM_020812.4</a>
UniProt ID:	<a href="#">Q96HP0</a>
Synonyms:	AOS2; ZIR1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)