

## Product datasheet for **HP213856**

### ZNF471 Human qPCR Primer Pair (NM\_020813)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57573
Forward Sequence:	GTTGGAGAACTATCAGAGCCTGG
Reverse Sequence:	GGCTTCTTGTCATCTCACTCGTC
ACCN:	<a href="#">NM_020813</a> , <a href="#">NM_020813.1</a> , <a href="#">NM_020813.2</a> , <a href="#">NM_020813.3</a> , <a href="#">BC125222</a> , <a href="#">BC016830</a> , <a href="#">BC028227</a> , <a href="#">BC125221</a> , <a href="#">BM991696</a> , <a href="#">NM_020813.4</a>
UniProt ID:	<a href="#">Q9BX82</a>
Synonyms:	ERP1; Z1971; Zfp78
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)