

## Product datasheet for **HP213896**

### SH3MD2 (SH3RF1) Human qPCR Primer Pair (NM\_020870)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57630
Forward Sequence:	GCTCACATCCAGACAAGTCCTC
Reverse Sequence:	TACCTGGATGGGTGCTCACTGCT
ACCN:	<a href="#">BC053671</a> , <a href="#">NM_020870</a> , <a href="#">NM_020870.1</a> , <a href="#">NM_020870.2</a> , <a href="#">NM_020870.3</a> , <a href="#">BC053671.1</a> , <a href="#">BC033203</a> , <a href="#">BC041023</a> , <a href="#">NM_020870.4</a>
UniProt ID:	<a href="#">Q7Z6I0</a>
Synonyms:	POSH; RNF142; SH3MD2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)