

Product datasheet for **HP213946**

MEAK7 (TLDC1) Human qPCR Primer Pair (NM_020947)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57707
Forward Sequence:	ATGGACAGCAGACGATCCCGAA
Reverse Sequence:	TTGGCTCTGCTGTGCCTTTCC
ACCN:	NM_020947 , NM_020947.1 , NM_020947.2 , NM_020947.3 , BC060844 , BC060844.1 , BC023251 , BG748167 , BI767086 , BM674034 , BM676714 , BM715013 , BM715579 , BM922004 , BQ009692 , BU955954 , NM_020947.4
UniProt ID:	Q6P9B6
Synonyms:	KIAA1609; mEAK-7; TLDC1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)