

## Product datasheet for **HP213981**

### MST1 Human qPCR Primer Pair (NM\_020998)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4485
Forward Sequence:	TGGTGCTACACGATGGACCCAA
Reverse Sequence:	GCCACACTTCTCAAAGTGCACC
ACCN:	<a href="#">NM_020998</a> , <a href="#">NM_020998.1</a> , <a href="#">NM_020998.2</a> , <a href="#">NM_020998.3</a> , <a href="#">BC048330</a> , <a href="#">BC048330.1</a> , <a href="#">BC044862</a> , <a href="#">BM672747</a>
UniProt ID:	<a href="#">P26927</a>
Synonyms:	D3F15S2; DNF15S2; HGFL; MSP; NF15S2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)