

Product datasheet for **HP213988**

LD78 beta (CCL3L1) Human qPCR Primer Pair (NM_021006)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6349
Forward Sequence:	ACTTTGAGACGAGCAGCCAGTG
Reverse Sequence:	TTTCTGGACCCACTCCTCACTG
ACCN:	<u>BC027888</u> , <u>NM_021006</u> , <u>NM_021006.1</u> , <u>NM_021006.3</u> , <u>NM_021006.4</u> , <u>NM_021006.5</u> , <u>BC027888.1</u> , <u>BC107710</u> , <u>BC146914</u>
UniProt ID:	<u>P16619</u>
Synonyms:	464.2; D17S1718; G0S19-2; LD78; LD78-beta(1-70); LD78BETA; MIP1AP; SCYA3L; SCYA3L1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)