

Product datasheet for **HP213997**

H3FI (HIST1H3F) Human qPCR Primer Pair (NM_021018)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8968
Forward Sequence:	TCGCAAGCTACCATTCCAGCGT
Reverse Sequence:	CAGGTTGGTGTCCTCAAAGAGC
ACCN:	<u>NM_021018</u> , <u>NM_021018.1</u> , <u>NM_021018.2</u> , <u>BC067492</u> , <u>BC062305</u> , <u>BC069818</u> , <u>BC096131</u> , <u>NM_021018.3</u>
UniProt ID:	<u>P68431</u>
Synonyms:	H3/i; H3C1; H3C2; H3C3; H3C4; H3C6; H3C8; H3C10; H3C11; H3C12; H3FI; HIST1H3F
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)