

## Product datasheet for **HP214006**

### RPL36A Human qPCR Primer Pair (NM\_021029)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6173
Forward Sequence:	AGGGCAAGGATTCTCTGTACGC
Reverse Sequence:	CAACGCACTCAAGCCTTAGCAC
ACCN:	<u><a href="#">NM_021029</a></u> , <u><a href="#">NM_021029.1</a></u> , <u><a href="#">NM_021029.2</a></u> , <u><a href="#">NM_021029.3</a></u> , <u><a href="#">NM_021029.4</a></u> , <u><a href="#">NM_021029.5</a></u> , <u><a href="#">BC062219</a></u> , <u><a href="#">BC001781</a></u> , <u><a href="#">BC031015</a></u> , <u><a href="#">BC070204</a></u> , <u><a href="#">BF211310</a></u> , <u><a href="#">BG504049</a></u> , <u><a href="#">BG943625</a></u> , <u><a href="#">BP245950</a></u> , <u><a href="#">NM_021029.6</a></u>
UniProt ID:	<u><a href="#">P83881</a></u>
Synonyms:	L36A; L44L; MIG6; RPL44
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)