

Product datasheet for **HP214009**

IFITM3 Human qPCR Primer Pair (NM_021034)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10410
Forward Sequence:	CTGGGCTTCATAGCATTGCGCT
Reverse Sequence:	AGATGTTCAAGCACTTGGCGGT
ACCN:	NM_021034 , NM_021034.1 , NM_021034.2 , BC006794 , BC008417 , BC022439 , BC070243 , BQ937265 , BT006892 , NM_021034.3
UniProt ID:	Q01628
Synonyms:	1-8U; DSPA2b; IP15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)