

Product datasheet for **HP214021**

H3C14 Human qPCR Primer Pair (NM_021059)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	126961
Forward Sequence:	GCGCTACCAGAAGTCCACGGA
Reverse Sequence:	AGGTCCGTCTTAAAGTCCTGCG
ACCN:	<u>NM_021059</u> , <u>NM_021059.1</u> , <u>NM_021059.2</u> , <u>BC074969</u> , <u>BC015544</u> , <u>BC130635</u> , <u>BC130637</u> , <u>BC148478</u> , <u>BC153074</u> , <u>NM_021059.3</u>
UniProt ID:	<u>Q71DI3</u>
Synonyms:	H3; H3.2; H3/M; H3C13; H3C15; H3F2; H3FM; H3FN; HIST2H3C
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)