

Product datasheet for **HP214024**

HIST1H2AK (HIST1H2AG) Human qPCR Primer Pair (NM_021064)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8969
Forward Sequence:	GCGACAACAAGAAGACCCGCAT
Reverse Sequence:	TGTGCGATGGTGACTTTGCCCA
ACCN:	<u>NM_021064</u> , <u>NM_021064.1</u> , <u>NM_021064.2</u> , <u>NM_021064.3</u> , <u>NM_021064.4</u> , <u>BC016677</u> , <u>BC016677.1</u> , <u>BC067782</u> , <u>NM_021064.5</u>
UniProt ID:	<u>P0C0S8</u>
Synonyms:	H2A.1b; H2A/p; H2AC13; H2AC15; H2AC16; H2AC17; H2AFP; H2AG; HIST1H2AG; pH2A/f
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)