

## Product datasheet for **HP214028**

### IFNA4 Human qPCR Primer Pair (NM\_021068)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3441
Forward Sequence:	GTTCCAGAAGGCTCAAGCCATC
Reverse Sequence:	TAGGAGGCTCTGTCCCAAGCA
ACCN:	<u><a href="#">NM_021068</a></u> , <u><a href="#">NM_021068.1</a></u> , <u><a href="#">NM_021068.2</a></u> , <u><a href="#">BC113640</a></u> , <u><a href="#">BC074965</a></u> , <u><a href="#">BC074966</a></u> , <u><a href="#">BC113642</a></u> , <u><a href="#">NM_021068.3</a></u>
UniProt ID:	<u><a href="#">P05014</a></u>
Synonyms:	IFN-alpha4a; INFA4
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)