

Product datasheet for **HP214031**

HCN1 Human qPCR Primer Pair (NM_021072)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	348980
Forward Sequence:	ACGAGAAGGAGCCGTGGGTAAA
Reverse Sequence:	ACGTCCTTTGGTCAGCAGGCAA
ACCN:	<u>BD139755</u> , <u>NM_021072</u> , <u>NM_021072.1</u> , <u>NM_021072.2</u> , <u>NM_021072.3</u> , <u>BC156454</u> , <u>BC172530</u> , <u>BM931647</u> , <u>BX537747</u> , <u>NM_021072.4</u>
UniProt ID:	<u>O60741</u>
Synonyms:	BCNG-1; BCNG1; DEE24; EIEE24; GEFSP10; HAC-2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)