

## Product datasheet for **HP214037**

### NMT1 Human qPCR Primer Pair (NM\_021079)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	4836
Forward Sequence:	GTTGCTCCAGTTCTGATCCGAG
Reverse Sequence:	TTCCGTGGGTTTAGGGACCGAT
ACCN:	<a href="#">NM_021079</a> , <a href="#">NM_021079.1</a> , <a href="#">NM_021079.2</a> , <a href="#">NM_021079.3</a> , <a href="#">NM_021079.4</a> , <a href="#">BC006538</a> , <a href="#">BC006538.2</a> , <a href="#">BC006569</a> , <a href="#">BC007258</a> , <a href="#">BC008312</a> , <a href="#">BC008579</a> , <a href="#">BC109226</a> , <a href="#">BU849351</a> , <a href="#">NM_021079.5</a>
UniProt ID:	<a href="#">P30419</a>
Synonyms:	NMT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH <sub>2</sub> O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)