

Product datasheet for **HP214043**

MTMR3 Human qPCR Primer Pair (NM_021090)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8897
Forward Sequence:	GTGCAACAACGCCAAGGAGAGA
Reverse Sequence:	CAGCATCAGGTTACGCACATGG
ACCN:	<u>NM_021090</u> , <u>NM_021090.1</u> , <u>NM_021090.2</u> , <u>NM_021090.3</u> , <u>BC141654</u> , <u>BC141656</u> , <u>BC142713</u> , <u>BC148216</u> , <u>BC151217</u> , <u>BC152455</u> , <u>BI462128</u> , <u>BU733010</u> , <u>NM_021090.4</u>
UniProt ID:	<u>Q13615</u>
Synonyms:	FYVE-DSP1; ZFYVE10
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)