

Product datasheet for **HP214049**

SPINT2 Human qPCR Primer Pair (NM_021102)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10653
Forward Sequence:	GTGCCTCAAGAAATGTGCCACTG
Reverse Sequence:	GGAGTGGTCTTCAGAATCCTGC
ACCN:	NM_021102 , NM_021102.1 , NM_021102.2 , NM_021102.3 , BC012868 , BC001668 , BC007705 , BC011951 , BC011955 , BT020115 , BT020116 , BT020117 , NM_021102.4
UniProt ID:	O43291
Synonyms:	DIAR3; HAI-2; HAI2; Kop; PB
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)