

Product datasheet for **HP214054**

COL14A1 Human qPCR Primer Pair (NM_021110)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7373
Forward Sequence:	CACAAACCTCCTCAGCGGAATG
Reverse Sequence:	GGCTTGGAGATTGGTAACACCC
ACCN:	<u>BC014640</u> , <u>NM_021110</u> , <u>NM_021110.1</u> , <u>NM_021110.2</u> , <u>NM_021110.3</u> , <u>BC014640.1</u> , <u>BC036192</u> , <u>BC083495</u> , <u>BC140893</u> , <u>BQ717750</u> , <u>BX472674</u> , <u>NM_021110.4</u>
UniProt ID:	<u>Q05707</u>
Synonyms:	UND
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)