

Product datasheet for **HP214058**

Adenylate cyclase 1 (ADCY1) Human qPCR Primer Pair (NM_021116)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	107
Forward Sequence:	TCGCGTCAACAGGTACATCAGC
Reverse Sequence:	CGCTGGTAAATACTCGTCCTG
ACCN:	BC036858 , NM_021116 , NM_021116.1 , NM_021116.2 , BC041473 , BC169213 , BC169214 , BU181704 , NM_021116.4
UniProt ID:	Q08828
Synonyms:	AC1; DFNB44
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)