

Product datasheet for **HP214060**

Cylicin 1 (CYLC1) Human qPCR Primer Pair (NM_021118)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	1538
Forward Sequence:	CAGTGAATCAAGCAGAAAATCATGG
Reverse Sequence:	GAGCTGGTTTCTGGCCTTCTTC
ACCN:	<u>NM_021118</u> , <u>NM_021118.1</u> , <u>NM_021118.2</u> , <u>BC126461</u> , <u>BC021669</u> , <u>BC038787</u> , <u>BC058019</u> , <u>BC070355</u> , <u>BX110674</u> , <u>NM_021118.3</u>
UniProt ID:	<u>P35663</u>
Synonyms:	CYCL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)