

Product datasheet for **HP214081**

COTLI Human qPCR Primer Pair (NM_021149)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	23406
Forward Sequence:	AAATATGACGGCTCCACCATCGT
Reverse Sequence:	TGGACCTCTTGCTCATGGCATC
ACCN:	<u>NM_021149</u> , <u>NM_021149.1</u> , <u>NM_021149.2</u> , <u>BC010039</u> , <u>BC010039.2</u> , <u>BC010884</u> , <u>BC016702</u> , <u>BC020445</u> , <u>BC042970</u> , <u>BC053682</u> , <u>BT006968</u> , <u>BX395081</u> , <u>NM_021149.5</u>
UniProt ID:	<u>Q14019</u>
Synonyms:	CLP
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.

