

## Product datasheet for **HP214083**

### CROT Human qPCR Primer Pair (NM\_021151)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	54677
Forward Sequence:	CTAGTGAGGAGCGAACTCGATG
Reverse Sequence:	CCTCTGGTGTTACATGTGGACTG
ACCN:	<u><a href="#">NM_021151</a></u> , <u><a href="#">NM_021151.1</a></u> , <u><a href="#">NM_021151.2</a></u> , <u><a href="#">NM_021151.3</a></u> , <u><a href="#">BC039004</a></u> , <u><a href="#">BC039004.1</a></u> , <u><a href="#">BC006015</a></u> , <u><a href="#">BC012483</a></u> , <u><a href="#">BC019226</a></u> , <u><a href="#">BC051874</a></u> , <u><a href="#">NM_021151.4</a></u>
UniProt ID:	<u><a href="#">Q9UKG9</a></u>
Synonyms:	COT
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)