

## Product datasheet for **HP214112**

### CADM3 Human qPCR Primer Pair (NM\_021189)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57863
Forward Sequence:	CTCGGTGACATTCCAGGTTACC
Reverse Sequence:	GCCTAATCATCGCAGTTGGTGTG
ACCN:	<u><a href="#">NM_021189</a></u> , <u><a href="#">NM_021189.1</a></u> , <u><a href="#">NM_021189.2</a></u> , <u><a href="#">NM_021189.3</a></u> , <u><a href="#">NM_021189.4</a></u> , <u><a href="#">BC033819</a></u>
UniProt ID:	<u><a href="#">Q8N126</a></u>
Synonyms:	BlgR; IGSF4B; Necl-1; NECL1; synCAM3; TSLL1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)