

Product datasheet for **HP214130**

CARD12 (NLRC4) Human qPCR Primer Pair (NM_021209)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	58484
Forward Sequence:	AGGTCCCACAACCTCGTCAAGCT
Reverse Sequence:	TGCTCACACGATTCCCGCCAA
ACCN:	<u>NM_021209</u> , <u>NM_021209.1</u> , <u>NM_021209.2</u> , <u>NM_021209.3</u> , <u>NM_021209.4</u> , <u>BC031555</u> , <u>BC031555.1</u>
UniProt ID:	<u>Q9NPP4</u>
Synonyms:	AIFEC; CARD12; CLAN; CLAN1; CLANA; CLANB; CLANC; CLAND; CLR2.1; FCAS4; IPAF
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)