

Product datasheet for **HP214152**

EPS15R (EPS15L1) Human qPCR Primer Pair (NM_021235)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	58513
Forward Sequence:	CAGAAGTGCCAGGATGAGACTC
Reverse Sequence:	CTCGGTTTCAGATCGTCTTCTCTG
ACCN:	NM_021235 , NM_021235.1 , NM_021235.2 , BC028243 , BC037558 , BC051873 , BC067114 , BC131590 , BC142662 , BC142716 , BU075732
UniProt ID:	Q9UBC2
Synonyms:	EPS15R
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)