

Product datasheet for **HP214164**

LILRA5 Human qPCR Primer Pair (NM_021250)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	353514
Forward Sequence:	GCAGGTATGGTCAGAACCCAGT
Reverse Sequence:	GCGGATGAGATTCTCTACTGCG
ACCN:	<u>NM_021250</u> , <u>NM_021250.1</u> , <u>NM_021250.2</u> , <u>NM_021250.3</u> , <u>BC148658</u> , <u>BC153203</u> , <u>NM_021250.4</u>
UniProt ID:	<u>A6NI73</u>
Synonyms:	CD85; CD85F; ILT-11; ILT11; LILRB7; LIR-9; LIR9
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)