

Product datasheet for **HP214179**

KCNN2 Human qPCR Primer Pair (NM_021614)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	3781
Forward Sequence:	GCCTTATCAGTCTCTCCACGATC
Reverse Sequence:	CCAGTCATCTGCTCCATTGTCC
ACCN:	<u>NM_021614</u> , <u>NM_021614.1</u> , <u>NM_021614.2</u> , <u>NM_021614.3</u> , <u>BC117456</u> , <u>BC015371</u> , <u>BC040722</u> , <u>BC117454</u> , <u>BG477734</u> , <u>BG769522</u>
UniProt ID:	<u>Q9H2S1</u>
Synonyms:	hSK2; KCa2.2; SK2; SKCA2; SKCa 2
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH2O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)