

Product datasheet for **HP214190**

PDLIM2 Human qPCR Primer Pair (NM_021630)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	64236
Forward Sequence:	GAGAAGTGCAGTACCAGCATCG
Reverse Sequence:	GCATCTTCAGTTAGCCACACA
ACCN:	NM_021630 , NM_021630.2 , NM_021630.3 , NM_021630.4 , NM_021630.5 , BC071774 , BC071774.1 , BC021556 , BM674373 , BU624103 , BU737724
UniProt ID:	Q96JY6
Synonyms:	MYSTIQUE; SLIM
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)