

Product datasheet for **HP214194**

LGR7 (RXFP1) Human qPCR Primer Pair (NM_021634)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	59350
Forward Sequence:	CACTGTAAACGGTGTGGACGACT
Reverse Sequence:	TCAGGTGTTTCTGCCTCAAAGG
ACCN:	<u>NM_021634</u> , <u>NM_021634.1</u> , <u>NM_021634.2</u> , <u>NM_021634.3</u> , <u>BC112142</u> , <u>BC112142.1</u> , <u>BC113617</u> , <u>BC143556</u> , <u>BC143557</u> , <u>BC143558</u> , <u>BC143559</u> , <u>BX647975</u> , <u>BX647980</u> , <u>BX647985</u> , <u>NM_021634.4</u>
UniProt ID:	<u>Q9HBX9</u>
Synonyms:	LGR7; RXFP1
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	Store at -20 °C.
Stability:	The primer mix is stable for one year from date of shipping.

[View online »](#)