

Product datasheet for **HP214213**

SLC25A19 Human qPCR Primer Pair (NM_021734)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	60386
Forward Sequence:	GCCATACCAGCCGAAGGAAAGA
Reverse Sequence:	CTCCAACCTGTAGCCGCTTCTT
ACCN:	NM_021734 , NM_021734.1 , NM_021734.2 , NM_021734.3 , NM_021734.4 , BC001075 , BC001075.2 , BC005120 , BM455504 , BU177306 , BU607955
UniProt ID:	Q9HC21
Synonyms:	DNC; MCPHA; MUP1; THMD3; THMD4; TPC
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)