

## Product datasheet for **HP214218**

### PLAC1 Human qPCR Primer Pair (NM\_021796)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	10761
Forward Sequence:	GCTCCATGAGAGTAGCCAGCAA
Reverse Sequence:	CAAGGTGGACAATCGCAGTTGG
ACCN:	<u><a href="#">NM_021796</a></u> , <u><a href="#">NM_021796.1</a></u> , <u><a href="#">NM_021796.2</a></u> , <u><a href="#">NM_021796.3</a></u> , <u><a href="#">BC022335</a></u> , <u><a href="#">BC022335.1</a></u> , <u><a href="#">BC066327</a></u> , <u><a href="#">BG436529</a></u> , <u><a href="#">NM_021796.4</a></u>
UniProt ID:	<u><a href="#">Q9HBJ0</a></u>
Synonyms:	CT92; OOSP2B; OOSP2L
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)