

Product datasheet for **HP214230**

ELOVL5 Human qPCR Primer Pair (NM_021814)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	60481
Forward Sequence:	ACGTCTACCACCATGCCTCGAT
Reverse Sequence:	TGGAAGGGACTGACGACAAACC
ACCN:	NM_021814 , NM_021814.1 , NM_021814.2 , NM_021814.3 , NM_021814.4 , BC067123 , BC067123.1 , BC009838 , BC017270 , BC074503 , NM_021814.5
UniProt ID:	Q9NYP7
Synonyms:	dj483K16.1; HELO1; SCA38
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)