

Product datasheet for **HP214239**

CCDC90B Human qPCR Primer Pair (NM_021825)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	60492
Forward Sequence:	TCAAATGTCAGCCTGGATACTATC
Reverse Sequence:	TCCTGATAGCATCCAATGAGCC
ACCN:	<u>NM_021825</u> , <u>NM_021825.1</u> , <u>NM_021825.2</u> , <u>NM_021825.3</u> , <u>NM_021825.4</u> , <u>BC014573</u> , <u>BC017771</u> , <u>BC020783</u> , <u>BC032701</u> , <u>BC048795</u> , <u>BC067764</u> , <u>BC107755</u> , <u>BF131344</u> , <u>BG700162</u>
UniProt ID:	<u>Q9GZT6</u>
Synonyms:	MDS011; MDS025
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)