

Product datasheet for **HP214286**

H4C11 Human qPCR Primer Pair (NM_021968)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8363
Forward Sequence:	TGCGCGACAATATCCAGGGCAT
Reverse Sequence:	GTTCTCCAGGAACACCTTCAGC
ACCN:	<u>BC017361</u> , <u>NM_021968</u> , <u>NM_021968.1</u> , <u>NM_021968.2</u> , <u>NM_021968.3</u> , <u>BC017361.2</u> , <u>BC140249</u> , <u>BC146523</u> , <u>BC148271</u> , <u>NM_021968.4</u>
UniProt ID:	<u>P62805</u>
Synonyms:	dj160A22.2; H4-16; H4/e; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4C12; H4C13; H4C14; H4C15; H4F2iv; H4FE; HIST1H4J
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)