

Product datasheet for **HP214288**

MAP2K1IP1 (LAMTOR3) Human qPCR Primer Pair (NM_021970)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	8649
Forward Sequence:	GTGGTGAGTTTCATAGCCAGCAG
Reverse Sequence:	AGAAACTTCCACAACTTGCTCTCAG
ACCN:	<u>NM_021970</u> , <u>NM_021970.1</u> , <u>NM_021970.3</u> , <u>BC026245</u> , <u>BC026245.1</u> , <u>BC005340</u> , <u>BC110402</u> , <u>BI562630</u> , <u>NM_021970.4</u>
UniProt ID:	<u>Q9UHA4</u>
Synonyms:	MAP2K1IP1; MAPBP; MAPKSP1; MP1; PRO0633; Ragulator3
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)