

Product datasheet for **HP214291**

POLR2F Human qPCR Primer Pair (NM_021974)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	5435
Forward Sequence:	AGCATGCACCTATTCCAGTGGC
Reverse Sequence:	ACAGGGAGATGGAGGGAAGGAT
ACCN:	<u>NM_021974</u> , <u>NM_021974.1</u> , <u>NM_021974.2</u> , <u>NM_021974.3</u> , <u>NM_021974.4</u> , <u>BC003582</u> , <u>BC003582.1</u> , <u>BP195499</u> , <u>BU160667</u> , <u>BU634038</u> , <u>NM_021974.5</u>
UniProt ID:	<u>P61218</u>
Synonyms:	HRBP14.4; POLRF; RPABC2; RPABC14.4; RPB6; RPB14.4; RPC15
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)