

Product datasheet for **HP214295**

ST14 Human qPCR Primer Pair (NM_021978)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	6768
Forward Sequence:	AGGACTACGTGGAGATCAACGG
Reverse Sequence:	TCACTGGAGTCGTAGGAGAGGT
ACCN:	<u>NM_021978</u> , <u>NM_021978.1</u> , <u>NM_021978.2</u> , <u>NM_021978.3</u> , <u>BC030532</u> , <u>BC030532.1</u> , <u>BC005826</u> , <u>BC018146</u> , <u>BM791311</u> , <u>BQ059215</u> , <u>NM_021978.4</u>
UniProt ID:	<u>Q9Y5Y6</u>
Synonyms:	ARCI11; CAP3; HAI; MT-SP1; MTSP1; PRSS14; SNC19; TADG15; TMPRSS14
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)