

Product datasheet for **HP214302**

ZNF711 Human qPCR Primer Pair (NM_021998)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	7552
Forward Sequence:	GCGGAGGCTGAAGATGATGTTG
Reverse Sequence:	CATTCGGCTCTGGTCAAGCACT
ACCN:	<u>NM_021998</u> , <u>NM_021998.1</u> , <u>NM_021998.2</u> , <u>NM_021998.3</u> , <u>NM_021998.4</u> , <u>BC067294</u> , <u>BC067294.1</u> , <u>BC006349</u> , <u>BX648117</u>
UniProt ID:	<u>Q9Y462</u>
Synonyms:	CMPX1; dj75N13.1; MRX65; MRX97; Zfp711; ZNF4; ZNF5; ZNF6
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)