

## Product datasheet for **HP214322**

### SLC4A10 Human qPCR Primer Pair (NM\_022058)

#### Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	57282
Forward Sequence:	TCAGTCTGCTCCAGCCTGTGTT
Reverse Sequence:	GATGCTTCAGCACCTGGAGGAA
ACCN:	<u><a href="#">NM_022058</a></u> , <u><a href="#">NM_022058.1</a></u> , <u><a href="#">NM_022058.2</a></u> , <u><a href="#">NM_022058.3</a></u> , <u><a href="#">BC036415</a></u> , <u><a href="#">BC038298</a></u> , <u><a href="#">BC051190</a></u> , <u><a href="#">BC136269</a></u> , <u><a href="#">BC143714</a></u> , <u><a href="#">NM_022058.4</a></u>
UniProt ID:	<u><a href="#">Q6U841</a></u>
Synonyms:	NBCn2; NCBE
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)