

Product datasheet for **HP214326**

C10orf84 (FAM204A) Human qPCR Primer Pair (NM_022063)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	63877
Forward Sequence:	CCGTCCTCAAATGAAACCCAGTG
Reverse Sequence:	ACTCCTCCACAGCCTGGTCTAT
ACCN:	<u>NM_022063</u> , <u>NM_022063.1</u> , <u>NM_022063.2</u> , <u>BC023577</u> , <u>BM837313</u> , <u>BX538256</u> , <u>NM_022063.3</u>
UniProt ID:	<u>Q9H8W3</u>
Synonyms:	bA319I23.1; C10orf84
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)