

Product datasheet for **HP214328**

E2 230K (UBE2O) Human qPCR Primer Pair (NM_022066)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	63893
Forward Sequence:	AACATCCGCTCCAACGACCTCT
Reverse Sequence:	CCAGACTGTACCACACCGTAGA
ACCN:	<u>NM_022066</u> , <u>NM_022066.1</u> , <u>NM_022066.2</u> , <u>NM_022066.3</u> , <u>BC004525</u> , <u>BC007822</u> , <u>BC022237</u> , <u>BC025977</u> , <u>BC036820</u> , <u>BC051868</u> , <u>NM_022066.4</u>
UniProt ID:	<u>Q9C0C9</u>
Synonyms:	E2-230K
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)