

Product datasheet for **HP214330**

FAM38B (PIEZO2) Human qPCR Primer Pair (NM_022068)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	63895
Forward Sequence:	GACGGACACAACCTTTGAGCCTG
Reverse Sequence:	CTGGCTTTGTTGGGCACTCATTG
ACCN:	<u>NM_022068</u> , <u>NM_022068.2</u> , <u>NM_022068.1</u> , <u>BC148379</u> , <u>BF509086</u> , <u>BI829775</u> , <u>BM974712</u>
UniProt ID:	<u>Q9H5I5</u>
Synonyms:	C18orf30; C18orf58; DA3; DA5; DAIPT; FAM38B; FAM38B2; HsT748; HsT771; MWKS
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)