

Product datasheet for **HP214364**

Claspin (CLSPN) Human qPCR Primer Pair (NM_022111)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	63967
Forward Sequence:	TGAAAACCTCAGGCAAGCAGCCTA
Reverse Sequence:	GGCTGATAGGATGGAATCGTGG
ACCN:	<u>NM_022111</u> , <u>NM_022111.1</u> , <u>NM_022111.2</u> , <u>BC115025</u> , <u>BC018670</u> , <u>BC038991</u> , <u>BC062215</u> , <u>BC113116</u> , <u>BC115023</u> , <u>BC115024</u> , <u>BC115026</u> , <u>BC115027</u> , <u>BC137279</u> , <u>BC140789</u> , <u>BG109920</u> , <u>NM_022111.4</u>
UniProt ID:	<u>Q9HAW4</u>
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.



[View online »](#)