

Product datasheet for **HP214366**

PRDM16 Human qPCR Primer Pair (NM_022114)

Product data:

Product Type:	qPCR Primer Pairs
Gene ID:	63976
Forward Sequence:	CAGCCAATCTCACCAGACACCT
Reverse Sequence:	GTGGCACTTGAAAGGCTTCTCC
ACCN:	<u>NM_022114</u> , <u>NM_022114.1</u> , <u>NM_022114.2</u> , <u>NM_022114.3</u> , <u>BC110593</u> , <u>BC161614</u> , <u>BC168363</u> , <u>BQ897896</u> , <u>BX092842</u> , <u>NM_022114.4</u>
UniProt ID:	<u>Q9HAZ2</u>
Synonyms:	CMD1LL; KMT8F; LVNC8; MEL1; PFM13
Components:	1 vial of lyophilized qSTAR qPCR primer mix (1 nmol each primer, sufficient for 200 reactions). Before use, reconstitute the primer mix in 200 µl dH ₂ O to make a final concentration of 10 µM.
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	Store at -20°C.
Stability:	The primer mix is stable for one year from date of shipping.


[View online »](#)